

```

#!/usr/bin/env python
# -*- coding: utf-8 -*-
"""
This experiment was created using PsychoPy3 Experiment Builder (v3.2.4),
on Thu Feb 13 09:46:41 2020
If you publish work using this script the most relevant publication is:

    Peirce J, Gray JR, Simpson S, MacAskill M, Höchenberger R, Sogo H, Kastman E, Lindeløv JK. (2019)
    PsychoPy2: Experiments in behavior made easy Behav Res 51: 195.
    https://doi.org/10.3758/s13428-018-01193-y
"""

from __future__ import absolute_import, division

from psychopy import locale_setup
from psychopy import prefs
from psychopy import sound, gui, visual, core, data, event, logging, clock
from psychopy.constants import (NOT_STARTED, STARTED, PLAYING, PAUSED,
                                STOPPED, FINISHED, PRESSED, RELEASED, FOREVER)

import numpy as np # whole numpy lib is available, prepend 'np.'
from numpy import (sin, cos, tan, log, log10, pi, average,
                   sqrt, std, deg2rad, rad2deg, linspace, asarray)
from numpy.random import random, randint, normal, shuffle
import os # handy system and path functions
import sys # to get file system encoding

from psychopy.hardware import keyboard

# Ensure that relative paths start from the same directory as this script
_thisDir = os.path.dirname(os.path.abspath(__file__))
os.chdir(_thisDir)

# Store info about the experiment session
psychopyVersion = '3.2.4'
expName = 'demo' # from the Builder filename that created this script
expInfo = {'participant': '', 'session': '001'}
dlg = gui.DlgFromDict(dictionary=expInfo, sortKeys=False, title=expName)
if dlg.OK == False:
    core.quit() # user pressed cancel
expInfo['date'] = data.getDateStr() # add a simple timestamp
expInfo['expName'] = expName
expInfo['psychopyVersion'] = psychopyVersion

# Data file name stem = absolute path + name; later add .psyexp, .csv, .log, etc
filename = _thisDir + os.sep + u'data/%s_%s_%s' % (expInfo['participant'], expName, expInfo['date'])

# An ExperimentHandler isn't essential but helps with data saving
thisExp = data.ExperimentHandler(name=expName, version="",
                                extraInfo=expInfo, runtimeInfo=None,
                                originPath='/Users/yangzhihan/Desktop/projects/perception_lab/demo_program/demo.py',
                                savePickle=True, saveWideText=True,
                                dataFileName=filename)
# save a log file for detail verbose info
logFile = logging.LogFile(filename+'.log', level=logging.EXP)
logging.console.setLevel(logging.WARNING) # this outputs to the screen, not a file

endExpNow = False # flag for 'escape' or other condition => quit the exp
frameTolerance = 0.001 # how close to onset before 'same' frame

# Start Code - component code to be run before the window creation

# Setup the Window
win = visual.Window(

```

```

65 size=(1024, 768), fullscr=True, screen=0,
66 winType='pyglet', allowGUI=False, allowStencil=False,
67 monitor='testMonitor', color=[0,0,0], colorSpace='rgb',
68 blendMode='avg', useFBO=True,
69 units='height')
70 # store frame rate of monitor if we can measure it
71 explInfo['frameRate'] = win.getActualFrameRate()
72 if explInfo['frameRate'] != None:
73     frameDur = 1.0 / round(explInfo['frameRate'])
74 else:
75     frameDur = 1.0 / 60.0 # could not measure, so guess
76
77 # create a default keyboard (e.g. to check for escape)
78 defaultKeyboard = keyboard.Keyboard()
79
80 # Initialize components for Routine "instructions"
81 instructionsClock = core.Clock()
82 instructions_text = visual.TextStim(win=win, name='instructions_text',
83     text='Welcome!',
84     font='Arial',
85     pos=(0, 0), height=0.1, wrapWidth=None, ori=0,
86     color='white', colorSpace='rgb', opacity=1,
87     languageStyle='LTR',
88     depth=0.0);
89 instructions_keyresp = keyboard.Keyboard()
90
91 # Initialize components for Routine "show_condition"
92 show_conditionClock = core.Clock()
93 show_condition_text = visual.TextStim(win=win, name='show_condition_text',
94     text='default text',
95     font='Arial',
96     pos=(0, 0), height=0.1, wrapWidth=None, ori=0,
97     color='white', colorSpace='rgb', opacity=1,
98     languageStyle='LTR',
99     depth=0.0);
100
101 # Create some handy timers
102 globalClock = core.Clock() # to track the time since experiment started
103 routineTimer = core.CountdownTimer() # to track time remaining of each (non-slip) routine
104
105 # -----Prepare to start Routine "instructions"-----
106 # update component parameters for each repeat
107 instructions_keyresp.keys = []
108 instructions_keyresp.rt = []
109 # keep track of which components have finished
110 instructionsComponents = [instructions_text, instructions_keyresp]
111 for thisComponent in instructionsComponents:
112     thisComponent.tStart = None
113     thisComponent.tStop = None
114     thisComponent.tStartRefresh = None
115     thisComponent.tStopRefresh = None
116     if hasattr(thisComponent, 'status'):
117         thisComponent.status = NOT_STARTED
118 # reset timers
119 t = 0
120 _timeToFirstFrame = win.getFutureFlipTime(clock="now")
121 instructionsClock.reset(-_timeToFirstFrame) # t0 is time of first possible flip
122 frameN = -1
123 continueRoutine = True
124
125 # -----Run Routine "instructions"-----
126 while continueRoutine:
127     # get current time
128     t = instructionsClock.getTime()
129     tThisFlip = win.getFutureFlipTime(clock=instructionsClock)
130     tThisFlipGlobal = win.getFutureFlipTime(clock=None)
131     frameN = frameN + 1 # number of completed frames (so 0 is the first frame)

```

component

initialize
all
components

component 2

these attributes are added
not assigned

2 styles of timers
one using conventional timer
the other using screen refresh

```

132 # update/draw components on each frame
133
134 # *instructions_text* updates
135 if instructions_text.status == NOT_STARTED and tThisFlip >= 0.0-frameTolerance:
136     # keep track of start time/frame for later
137     instructions_text.frameNStart = frameN # exact frame index
138     instructions_text.tStart = t # local t and not account for scr refresh
139     instructions_text.tStartRefresh = tThisFlipGlobal # on global time
140     win.timeOnFlip(instructions_text, 'tStartRefresh') # time at next scr refresh
141     instructions_text.setAutoDraw(True)
142
143 # *instructions_keyresp* updates
144 waitOnFlip = False
145 if instructions_keyresp.status == NOT_STARTED and tThisFlip >= 0.0-frameTolerance:
146     # keep track of start time/frame for later
147     instructions_keyresp.frameNStart = frameN # exact frame index
148     instructions_keyresp.tStart = t # local t and not account for scr refresh
149     instructions_keyresp.tStartRefresh = tThisFlipGlobal # on global time
150     win.timeOnFlip(instructions_keyresp, 'tStartRefresh') # time at next scr refresh
151     instructions_keyresp.status = STARTED
152     # keyboard checking is just starting
153     waitOnFlip = True
154     win.callOnFlip(instructions_keyresp.clock.reset) # t=0 on next screen flip
155     win.callOnFlip(instructions_keyresp.clearEvents, eventType='keyboard') # clear events on next screen flip
156 if instructions_keyresp.status == STARTED and not waitOnFlip:
157     theseKeys = instructions_keyresp.getKeys(keyList=['space'], waitRelease=False)
158     if len(theseKeys):
159         theseKeys = theseKeys[0] # at least one key was pressed
160
161     # check for quit:
162     if "escape" == theseKeys:
163         endExpNow = True
164         instructions_keyresp.keys = theseKeys.name # just the last key pressed
165         instructions_keyresp.rt = theseKeys.rt
166         # a response ends the routine
167         continueRoutine = False
168
169 # check for quit (typically the Esc key)
170 if endExpNow or defaultKeyboard.getKeys(keyList=["escape"]):
171     core.quit()
172
173 # check if all components have finished
174 if not continueRoutine: # a component has requested a forced-end of Routine
175     break
176 continueRoutine = False # will revert to True if at least one component still running
177 for thisComponent in instructionsComponents:
178     if hasattr(thisComponent, "status") and thisComponent.status != FINISHED:
179         continueRoutine = True
180         break # at least one component has not yet finished
181
182 # refresh the screen
183 if continueRoutine: # don't flip if this routine is over or we'll get a blank screen
184     win.flip()
185
186 # -----Ending Routine "instructions"-----
187 for thisComponent in instructionsComponents:
188     if hasattr(thisComponent, "setAutoDraw"):
189         thisComponent.setAutoDraw(False)
190 thisExp.addData('instructions_text.started', instructions_text.tStartRefresh)
191 thisExp.addData('instructions_text.stopped', instructions_text.tStopRefresh)
192 # check responses
193 if instructions_keyresp.keys in [None]: # No response was made
194     instructions_keyresp.keys = None
195 thisExp.addData('instructions_keyresp.keys', instructions_keyresp.keys)
196 if instructions_keyresp.keys != None: # we had a response
197     thisExp.addData('instructions_keyresp.rt', instructions_keyresp.rt)

```

why set twice? timeOnFlip is more accurate

Determines whether the stimulus should be automatically drawn on every frame flip.

process keyboard inputs

decide whether this Routine should end

record data

```

199 thisExp.addData('instructions_keyresp.started', instructions_keyresp.tStartRefresh)
200 thisExp.addData('instructions_keyresp.stopped', instructions_keyresp.tStopRefresh)
201 thisExp.nextEntry()
202 # the Routine "instructions" was not non-slip safe, so reset the non-slip timer
203 routineTimer.reset()
204
205 # set up handler to look after randomisation of conditions etc
206 trials = data.TrialHandler(nReps=1, method='random',
207     extraInfo=expInfo, originPath=-1,
208     trialList=data.importConditions('csvs/noise_conditions.csv'),
209     seed=None, name='trials')
210 thisExp.addLoop(trials) # add the loop to the experiment
211 thisTrial = trials.trialList[0] # so we can initialise stimuli with some values
212 # abbreviate parameter names if possible (e.g. rgb = thisTrial.rgb)
213 if thisTrial != None:
214     for paramName in thisTrial:
215         exec('{} = thisTrial[paramName]'.format(paramName))
216
217 for thisTrial in trials: loop through all the trials
218     currentLoop = trials
219     # abbreviate parameter names if possible (e.g. rgb = thisTrial.rgb)
220     if thisTrial != None:
221         for paramName in thisTrial:
222             exec('{} = thisTrial[paramName]'.format(paramName))
223
224 # -----Prepare to start Routine "show_condition"-----
225 routineTimer.add(1.000000)
226 # update component parameters for each repeat
227 show_condition_text.setText(condition)
228 # keep track of which components have finished
229 show_conditionComponents = [show_condition_text]
230 for thisComponent in show_conditionComponents:
231     thisComponent.tStart = None
232     thisComponent.tStop = None
233     thisComponent.tStartRefresh = None
234     thisComponent.tStopRefresh = None
235     if hasattr(thisComponent, 'status'):
236         thisComponent.status = NOT_STARTED
237
238 # reset timers
239 t = 0
240 _timeToFirstFrame = win.getFutureFlipTime(clock="now")
241 show_conditionClock.reset(-_timeToFirstFrame) # t0 is time of first possible flip
242 frameN = -1
243 continueRoutine = True
244
245 # -----Run Routine "show_condition"-----
246 while continueRoutine and routineTimer.getTime() > 0: this loop is forced stopped by this
247     # get current time
248     t = show_conditionClock.getTime()
249     tThisFlip = win.getFutureFlipTime(clock=show_conditionClock)
250     tThisFlipGlobal = win.getFutureFlipTime(clock=None)
251     frameN = frameN + 1 # number of completed frames (so 0 is the first frame)
252     # update/draw components on each frame
253
254     # *show_condition_text* updates
255     if show_condition_text.status == NOT_STARTED and tThisFlip >= 0.0-frameTolerance:
256         # keep track of start time/frame for later
257         show_condition_text.frameNStart = frameN # exact frame index
258         show_condition_text.tStart = t # local t and not account for scr refresh
259         show_condition_text.tStartRefresh = tThisFlipGlobal # on global time
260         win.timeOnFlip(show_condition_text, 'tStartRefresh') # time at next scr refresh
261         show_condition_text.setAutoDraw(True)
262     if show_condition_text.status == STARTED:
263         # is it time to stop? (based on global clock, using actual start)
264         if tThisFlipGlobal > show_condition_text.tStartRefresh + 1-frameTolerance:
265             # keep track of stop time/frame for later

```

not used

```

266 show_condition_text.tStop = t # not accounting for scr refresh
267 show_condition_text.frameNStop = frameN # exact frame index
268 win.timeOnFlip(show_condition_text, 'tStopRefresh') # time at next scr refresh
269 show_condition_text.setAutoDraw(False)
270
271 # check for quit (typically the Esc key)
272 if endExpNow or defaultKeyboard.getKeys(keyList=["escape"]):
273     core.quit()
274
275 # check if all components have finished
276 if not continueRoutine: # a component has requested a forced-end of Routine
277     break
278 continueRoutine = False # will revert to True if at least one component still running
279 for thisComponent in show_conditionComponents:
280     if hasattr(thisComponent, "status") and thisComponent.status != FINISHED:
281         continueRoutine = True
282         break # at least one component has not yet finished
283
284 # refresh the screen
285 if continueRoutine: # don't flip if this routine is over or we'll get a blank screen
286     win.flip()
287
288 # -----Ending Routine "show_condition"-----
289 for thisComponent in show_conditionComponents:
290     if hasattr(thisComponent, "setAutoDraw"):
291         thisComponent.setAutoDraw(False)
292 trials.addData('show_condition_text.started', show_condition_text.tStartRefresh)
293 trials.addData('show_condition_text.stopped', show_condition_text.tStopRefresh)
294 thisExp.nextEntry()
295
296 # completed 1 repeats of 'trials'
297
298
299 # Flip one final time so any remaining win.callOnFlip()
300 # and win.timeOnFlip() tasks get executed before quitting
301 win.flip()
302
303 # these shouldn't be strictly necessary (should auto-save)
304 thisExp.saveAsWideText(filename+'.csv')
305 thisExp.saveAsPickle(filename)
306 logging.flush()
307 # make sure everything is closed down
308 thisExp.abort() # or data files will save again on exit
309 win.close()
310 core.quit()

```

update screen

stop all components
from drawing

record start and stop

proceed to next entry